

## MW150 MAX



pH / ORP / Temperature Laboratory Bench Meter

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

#### **MW150 MAX** pH / ORP / Temperature Laboratory Bench Meter

**MW150 MAX** is an advanced pH/ORP/Temp microprocessor-based bench meter. It is ideal for students and technicians who need fast and reliable measurements.

This meter is provided with a series of new diagnostic features which add an entirely new dimension to the measurement of pH, by allowing the user to dramatically improve the reliability of the measurement:

- Up to 3-point automatic calibration with 7 standard calibration buffers
- Automatic or manual temperature compensation
- Built-in rechargeable battery with 8 hours battery life
- Auto-off feature to preserve battery energy
- Battery charger with battery monitor
- Dedicated GLP key
- Alphanumeric LCD displayed messages for user friendly, intuitive information/warning/error messages
- Internal clock and date to keep track of different time-dependent functions (calibration, timestamp, calibration time out)
- Dedicated GLP key

**Accessories:**

- **MA9001** pH 1.68 buffer solution, 230 mL bottle
- **MA9004** pH 4.01 buffer solution, 230 mL bottle
- **MA9006** pH 6.86 buffer solution, 230 mL bottle
- **MA9007** pH 7.01 buffer solution, 230 mL bottle
- **MA9009** pH 9.18 buffer solution, 230 mL bottle
- **MA9010** pH 10.01 buffer solution, 230 mL bottle
- **MA9012** Refilling solution for pH electrode, 230 mL bottle
- **MA9015** Electrode storage solution, 230 mL bottle
- **MA9016** Electrode cleaning solution, 230 mL bottle
- **MA9112** pH 12.45 buffer solution, 230 mL bottle
- **MA9315** Electrode Holder
- **MA917B/1** Glass body, double junction refillable pH electrode
- **MA831R** Temperature probe

**Ordering information:**

MW150 is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
- MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution
- M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution
- M10016 Sachet Electrode Cleaning Solution
- 12 VDC Adapter
- Instruction manual