

MW106 MAX



pH / ORP / Temperature MAX Portable Meter

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

MW106 MAX pH / ORP / Temperature Professional Portable Meter

The Milwaukee **MW106 MAX** meter combines all the features of a benchtop unit into a portable, IP67 rated meter.

The instruments are provided with a series of new diagnostic features for improved pH measurement reliability:

- IP67 waterproof casing
- pH extended range: -2.00 to 20.00 pH/ -2.000 to 20.000 pH up to 5-point automatic pH calibration
- 7 standard calibration buffers (pH 1.68,4.01, 6.86, 7.01, 9.18, 10.01 and 12.45) and two custom buffers
- log space for up to 1000 records
- alphanumeric LCD displayed messages for user friendly, intuitive information / warning / error messages
- auto-off feature to preserve battery life
- internal clock and date to keep track of different
- time-dependent functions (calibration timestamp, calibration time out)
- dedicated GLP key

Accessories:

- **MA906BR/1** pH//Temp probe with 1 m cable
- **M10004B** pH 4.01 buffer solution 20 mL sachet (25 pcs)
- **M10006B** pH 6.86 buffer solution 20 mL sachet (25 pcs)
- **M10007B** pH 7.01 buffer solution 20 mL sachet (25 pcs)
- **M10010B** pH 10.01 buffer solution 20 mL sachet (25 pcs)
- **MA9004** pH 4.01 buffer solution, 230 mL bottle
- **MA9006** pH 6.86 buffer solution, 230 mL bottle
- **MA9007** pH 7.01 buffer solution, 230 mL bottle
- **MA9009** pH 9.18 buffer solution, 230 mL bottle
- **MA9010** pH 10.01 buffer solution, 230 mL bottle
- **MA9015** Electrode storage solution, 230 mL
- **MA9016** Electrode cleaning solution, 230 mL
- **M10000B** Electrode rinse solution, 20 mL (25 pcs)

Ordering information:

MW106 is supplied complete with MA906BR/1 pH/Temp probe, 20 mL pH 4.01 and 7.01 sachet of calibration solution, 2x20 mL sachet of electrode cleaning solutions, Micro USB cable, 9V battery, instructions, all in a rugged carrying case.