

## TH310

PRODUCT  
Image  
Unavailable

TH310



Pocket Thermometer with Stainless Steel Probe

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

#### **TH310 Pocket Thermometer with Stainless Steel Probe**

Scientists and laboratory technicians rely on the accuracy of their thermometers when performing routine measurements. For this reason, Milwaukee developed the TH310. This palm-sized unit is a highly accurate thermometer that is destined to make glass thermometers obsolete.

TH310's features:

- Easy-to-read LCD displays clear reading at any angle
- Wide measurements range: -50.0°C to +150.0°C
- Fast response with a high-tech NTC thermistor sensor
- Up to 3000 hours of battery life
- Stainless steel probe resists aggressive environments

## Thermometers: TH310

---

### Ordering information:

TH300 is supplied with stainless steel probe with 1 meter cable, batteries and instruction manual.