

MW600 PRO



PRO Portable Dissolved Oxygen Meter

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

MW600 PRO Portable Dissolved Oxygen Meter

The MW600 is a Portable Dissolved Oxygen meter ideal for use in school laboratories. Dissolved Oxygen measurements are also very important in fish farms and ponds, where Oxygen levels are continuously monitored to obtain optimal reproduction.

The MW600 calibrates easily in 2 points (at 100% saturated air and in 0 Oxygen solution) and has Automatic Temperature Compensation which guarantees the highest accuracy. The low battery warning, easy to replace screw on cap membranes make this meter very simple to operate. Rugged Carrying Case (Optional) provides handy on-site meter calibration and measurements.

MW600 is supplied complete with a MA840 D.O. polarographic probe with 3 m cable, calibration screwdriver, 2 spare membranes, MA7040 (20 mL) electrolyte solution, battery and instructions.

Accessories:

- **MA9070** Zero Oxygen Solution, 230 ml bottle
- **MA9071** Refilling Electrolyte Solution, 230 ml bottle
- **MA840** Polarographic D.O. probe with 3 meters cable
- **MA841** Spare membrane (5 pcs)

Ordering information:

MW600 is supplied complete with MA840 probe, 2 spare membranes, 20 mL bottle of electrolyte solution, calibration screwdriver, 9V battery and instructions.