

MW500 PRO

PRODUCT
Image
Unavailable



PRO Portable ORP Meter

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

MW500 PRO Portable ORP Meter

PRO portable meters with no frills! Milwaukee's PRO durable meters for quick and reliable measurements.

Milwaukee's PRO meters are manufactured to be easy to use, practical and accurate. Ideal for the classroom, laboratory or for general field use.

MW500 performs ORP measurements with a range of ± 1000 mV.

Accessories:

- **SE300** Platinum ORP electrode with 1 m cable
- **MA9015** Electrode storage solution, 230 ml bottle
- **MA9016** Electrode cleaning solution, 230 ml bottle
- **M10000B** Electrode rinse solution, 20 ml sachet (25 pcs)
- **MA9020** 200-275 mV ORP solution, 230 ml bottle

Ordering information:

MW500 is supplied complete with a SE300 platinum electrode, 9V battery and instructions.