

MW101 PRO

PRODUCT
Image
Unavailable



PRO Portable pH Meter with 0.01 pH resolution

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

MW101 PRO Portable pH Meter with 0.01 pH resolution

PRO portable meters with no frills! Milwaukee's durable meters for quick and reliable measurements.

Milwaukee's PRO meters are manufactured to be easy to use, practical and accurate. Ideal for the classroom, laboratory or for general field use.

MW101 performs pH measurements with a 0.01 pH resolution and with manual temperature compensation. The meter is supplied with pH electrode and calibration solutions.

Accessories:

- **SE220** pH electrode with BNC connector and 1 m cable
- **M10004B** pH 4.01 buffer solution 20 ml sachet (25 pcs)
- **M10007B** pH 7.01 buffer solution 20 ml sachet (25 pcs)
- **M10010B** pH 10.01 buffer solution 20 ml sachet (25 pcs)
- **MA9004** pH 4.01 buffer solution, 230 ml bottle
- **MA9007** pH 7.01 buffer solution, 230 ml bottle
- **MA9010** pH 10.01 buffer solution, 230 ml bottle
- **MA9015** Electrode storage solution, 230 ml bottle
- **MA9016** Electrode cleaning solution, 230 ml bottle
- **M10000B** Electrode rinse solution, 20 ml sachet (25 pcs)

Ordering information:

MW101 is supplied complete with a SE220 pH electrode, pH 7.01 20 mL sachet of calibration solution, calibration screwdriver, 9V battery and instructions.